

AmpliSens® PCR & Real-Time PCR DIAGNOSTIC KITS

Oncological diseases / 腫瘍疾患							Number of tests	Shelf life*
Type & package	Cat.No	Name	Description					
🕒 📄 ↺	TR-01(RG,iQ,Mx,ABI)-CE	AmpliSens® Leukosis-Quantum M-bcr-FRT	RNA extraction, Reverse transcription and Real-Time PCR Kits for use with RG, iQ, Mx and ABIPrism 7000. QUANTITATIVE. Determination of the <i>bcr-abl</i> chimeric gene (<i>M-bcr</i> variant) mRNA and <i>abl</i> gene mRNA concentration	CE		50/100	12 months	

HIV and HIV-associated infections / HIVとHIV関連の感染症							Number of tests	Shelf life*
HIV								
Type & package	Cat.No	Name	Description					
🕒 📄 ↺	TR-V0-G(RG,iQ)-CE	AmpliSens® DNA-HIV-FRT	DNA extraction and Real-Time PCR kits for use with RG and iQ	RVO		100	6 months	
🕒 📄 ↺	R-V0-M (RG,iQ,Mx,Dt)-CE	AmpliSens® HIV-Monitor-FRT	Reverse Transcription and Real-Time PCR kits for use with RG, iQ, Mx, and Dt QUANTITATIVE. Determination of HIV RNA concentration in clinical samples	RVO		76	6 months	
🕒 📄 ↺	R-V0-R(RG,iQ,Mx)-CE	AmpliSens® RNA-HIV-FRT	Real-Time PCR kit for use with RG, iQ, and Mx	RVO		76	6 months	

Identification of drug resistant mutations							Number of tests	Shelf life*
Type & package	Cat.No	Name	Description					
🕒 📄 ↺	R-02(RG,iQ)-CE	AmpliSens® Genoscreen HLA B*5701-FRT	Real time PCR kit for use with RG and iQ Kit for identification of 5701 allele of locus B in Human Major Histocompatibility Complex (HLA B*5701)	CE		110	12 months	

HIV associated infections							Number of tests	Shelf life*
Type & package	Cat.No	Name	Description					
🕒 📄 ↺	R-F2(RG,iQ,Mx)-CE	AmpliSens® Pneumocystis jirovecii (carinii)-FRT	Real-Time PCR kit for use with RG, iQ, and Mx	CE		55	6 months	

«MultiPrime» series (Multiplex pathogen detection)							Number of tests	Shelf life*
Type & package	Cat.No	Name	Description					
🕒 📄 ↺	R-V50-4x (RG,iQ,Mx,Dt)-CE	AmpliSens® HCV/HBV/HIV-FL	Reverse transcription and Real-Time PCR kits for use with RG, iQ, and Dt Differentiation of Hepatitis C (HCV), B (HBV) and HIV-I 4 optical channels are required	CE		100	9 months	

Legend: 🕒FRT; 📄FEP; 📄EPh; 📄-kit of reagents for PCR or Real-Time PCR;
 ⇓-reaction mix pre-dispensed into individual PCR tubes, under wax;
 ↺-all reagents are in separate vials (non aliquoted) + Taq-F polymerase (<hot start>)
 (*)Manufacturer shelf life - we guarantee not less than 70% of this shelf life.

Number of tests includes controls and calibrators (for quantitative tests)

Instruments	RG	Rotor-Gene Q (QIAGEN, Germany), Rotor-Gene 3000, 6000 (Corbett Research, Australia)
	iQ	iCycler iQ/iQ5 (Bio-Rad, USA)
	Mx	Mx3000P/3005P (Stratagene, USA)
	SC	SmartCycler (Cepheid, USA)
	Dt	DTprime (DT-96) (DNA-Technology, Russia)
CE	CE-marked kits, comply with EU directives 93/42/EEC and 98/79/EC for medical products and IVD	
RVO	Kits for Research Use Only	